



**Abbott**

# SCALE UP COVID-19 ANTIBODY TESTING

With a rapid lateral flow,  
fingerstick IgG/IgM assay

FOR EXTERNAL USE

Product not available in all countries. Not approved for sale in the U.S.





# Introducing the Panbio™ COVID-19 IgG/IgM Rapid Test Device

# Product Overview

- Lateral Flow Rapid Test Device
- Detects IgG/IgM antibodies caused by SARS-CoV-2
- Antibody test are utilized to determine if someone was previously infected
- Separate test lines for IgG and IgM
- 20µL fingerstick blood sample
- Time to result: 10-20 minutes
- Simple test procedure



# Kit Contents

## Materials Provided:

- 25 Test Devices individually pouched
- 25 Specimen droppers (for fingerstick whole blood only)
- 1 Buffer (3mL/vial)
- 1 Instructions for Use



## Materials Required but Not Provided:

- Specimen collection equipment and containers
- Micropipette
- Lancet (for fingerstick whole blood only)
- Cotton wool or gauze pad (for fingerstick whole blood only)
- Centrifuge
- Timer
- Biohazard waste containers for sharps and non-sharps

## Intended Use

The Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) is an *in vitro* diagnostic rapid test for the qualitative detection of IgG and IgM antibodies to SARS-CoV-2 in human serum, plasma, venous and fingerstick whole blood.

The Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) is for professional use only and is intended to be used as an aid in the diagnosis of SARS-CoV-2 infection. The product may be used in any laboratory and non-laboratory environment that meets the requirements specified in the Instructions for Use and local regulation.


The test provides preliminary test results. Negative results will not preclude SARS-CoV-2 infection and they cannot be used as the sole basis for treatment or other management decision.

The test is not intended to be used as a donor screening test for SARS-CoV-2.


# Fingerstick Test Procedure

Consult Instructions for Use for complete procedure


1 PUNCTURE FINGER WITH STERILE LANCET



2 ADD 20  $\mu$ L SAMPLE




3 ADD 2 DROPS OF BUFFER



4 READ RESULTS AT 10 MINUTES

DO NOT READ AFTER 20 MINUTES



# Performance

**Sensitivity: 95.8%** (95%CI: 85.7%~99.5%)

**Specificity: 94.0%** (95%CI: 83.5%~98.7%)

		PCR Positive	Clinical Negative	Total
Panbio™ COVID-19 IgG/IgM Rapid Test Device Result	IgG or IgM Positive	46	3	49
	IgG and IgM Negative	2	47	49
	Total	48	50	98


# Specifications

<b>Sample Types</b>	Fingerstick whole blood, venous whole blood, serum/plasma
<b>Sample Size</b>	-Fingerstick and venipuncture whole blood: 20µL -Serum/plasma: 10µL
<b>Test Time</b>	10 minutes, do not read after 20 minutes
<b>CE Mark</b>	Yes
<b>Sensitivity</b>	95.8%
<b>Specificity</b>	94.0%
<b>Storage</b>	2°–30°C
<b>Shelf Life</b>	12 months
<b>Qty per kit</b>	25 tests
<b>Control</b>	Procedural control included within the test
<b>Cat No</b>	ICO-T402



# Diagnosis Overview: IgG/IgM

# General Viral Response In Patient



NOTE: this is a typical viral patient response and is not specific to COVID-19

**In the case of COVID-19 (SARS-CoV-2 virus) the duration of the different stages of infection is still under investigation.** COVID-19 is estimated to have a mean incubation period of 6.4 days.<sup>1</sup>

- Many patients can be asymptomatic during the incubation phase

1. Lai CC, Shih TP, Ko WC, et.al. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. Int J Antimicrob Agents. 2020; Mar 55(3): 105924. doi: <https://doi.org/10.1016/j.ijantimicag.2020.105924>

# Growing Need for Rapid Point-of-care Tests

*“Development of rapid and accurate point-of-care tests which perform well in field settings are especially useful...*

*This would markedly improve early detection and isolation of infected patients and, by extension, identification of contacts.”<sup>1</sup>*

**-Report of the WHO-China Joint Mission on Coronavirus Disease 2019**

# Using IgG/IgM to Track Outbreak Expansion<sup>1</sup>

- IgG/IgM result is an indication of exposure to the virus
- Cannot be used solely for treatment or decision but can provide an indication of contact or potential infection
- Quickly identifying potential cases is key when tracking expansion of outbreak

# COVID-19 Patient Response in Literature

- **Different sources note 80% of COVID-19 cases seem to be asymptomatic or with mild symptoms.**
  - Molecular Diagnostic used for patients with symptoms, and in some countries, only very sick patients.
  - Asymptomatic patients can spread the virus without knowing.
- **IgM appears to show up earlier compared to other viral infections, approximately day 3 to 5 following infection.<sup>1</sup>**
  - IgM is an important biomarker, specifically for asymptomatic patients.
  - IgM has a short duration, estimating 2 to 4 weeks, but more data specific to COVID-19 is needed.
- **IgG also shows up earlier during COVID-19 infection.**
  - Testing for IgG/IgM can be a very effective diagnostic complement on asymptomatic patients.
  - Help identify potential asymptomatic patients quickly, increase contact case finding and reinforce quarantine need.

1. Diagnosis and Treatment Protocol for Novel Coronavirus Pneumonia, (Trial Version 7), March 3, 2020;  
[http://www.kankyokansen.org/uploads/uploads/files/jsipc/protocol\\_V7.pdf](http://www.kankyokansen.org/uploads/uploads/files/jsipc/protocol_V7.pdf)



**Lateral Flow Allows Rapid  
Scale Up of Testing Capacity**

## Scaling Up with the Panbio COVID-19 IgG/IgM Test

- Make testing more widely available
- Address urgent testing shortfall
- Decentralized testing can help to alleviate overcrowding in hospitals
- Reserve central facilities for seriously ill
- Useful in drive-through, home visit, point-of-care, or lab settings

# Empower frontline workers



## **Panbio™** COVID-19 IgG/IgM Rapid Test Device







**Abbott**

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model.

120006510-03 04/20